MOBOTIX NurseAssist Smart Sensor

powered by Kepler Vision Technologies

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MOBOTIX c71 Kepler NurseAssist Smart Sensor

inkl. Kepler NurseAssist App (Standard)
& Extended Demo (90 Tage)

MOBOTIX c71 Kepler NurseAssist Smart Sensor – Extended Lizenz

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Kepler NurseAssist Extended Plug-In Softwarelizenz

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MOBOTIX Nurse Call & Reporting Software für c71 Kepler NurseAssist Smart Sensoren

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Integration in externes Nurse Call /Alarm Management System

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NurseAssist Integrations Guide #1

MOBOTIX HUB & NurseAssist Dashboard Plug-In

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+

+

MOBOTIX c71 NurseAssist Smart Sensor

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
**MOBOTIX NurseAssist Privacy
Whitepaper**

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The screenshot displays the 'Kepler NurseAssist Settings' window for a MOBOTIX c71 camera. The interface is organized into several sections:

- notifications**:
 - Enable in-bed notifications**: Checked (checkbox).
 - In-bed/not-in-bed notification duration**: Set to 30 seconds.
- MOBOTIX HUB Analytic Event**:
 - Enable**: Checked (checkbox).
 - MOBOTIX HUB Host Address**: 10.250.250.63.
 - MOBOTIX HUB Port Number**: 9090.
 - Camera Name**: 10.32.211.95.
- Advanced**:
 - Ignore areas Polygon points**: 7 x 5.

At the bottom of the settings window, there are four buttons: 'Set', 'Factory', 'Restore', and 'Close'.

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Profiles & Options	Value	Explanation
IP Notify Type	Custom Configuration	Predefined Configuration: "MxCC Alarm" sends predefined network messages to the MxCC alarm list. Acknowledge Required prompts the MxCC user to confirm the message. If the alarm is not acknowledged within the specified acknowledge time, the camera triggers a transmission error. Select Custom Configuration to see the extended configuration.
Destination Address	192.168.10.20:8000	Destination Addresses: Receiver IP address and port. Separate IP address and port using a colon. Enter one address per line.
	Parallel send to all	Send Order: Send notification to one or more destinations. Sequential and parallel will send a notification to each destination address. Send to next on error will stop after the first successful notification or will try the next address if unsuccessful.
Data Protocol	Raw TCP/IP	Transfer Protocol: Transfer notification data using these protocol headers.
	--next image--	Separator for Raw TCP/IP: Enter the separator for splitting several notification parts in Raw TCP/IP mode. Enter a unique string.
		CGI Path: Absolute CGI path beginning with '/'. This parameter allows using variables .
	*****	HTTP Authentication: User name and password for HTTP authentication separated by colon. Example: admin:meinam
Data Type	Plain text	Notification Data: Select type of IP notification data.
	[{"name":"\${ID.NAM},"event":"\${FPRENO},"eventName":"VM2","ip":"\${ID.ETHERNET}"]	Message: Message to include in Plain text notification data. When using HTTP protocol this text is used for QUERY_STRING in GET request. This parameter allows using variables .

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General Settings	Value	Explanation
Action Group	IQM VM2 Falling	Name: The name is purely informational.
	Enabled	Arming: Controls this action group: Enabled: activate the group. Off: deactivate the group. SI: group armed by signal input. CS: group armed by custom signal as defined in General Event Settings
	(No time table)	Time Table: Time table for this action profile (Time Tables).
Event Selection	<ul style="list-style-type: none"> (Message: KeplerNurseAssist) (Message: Fall_Detected) (Message: Not_In_Bed) (Message: In_Bed) (Message: Situation_Normal) 	Event Selection: Select the events which will trigger the actions below. Use [Ctrl] Click to select more than one event. Events in parentheses need to be activated first.
Action Details	5	Action Deadline: Time to wait (0..3600 s) before a new action can take place.
	Simultaneously	Action Chaining: Choose how the status of each subaction influences the execution of all others. Simultaneously: All actions are executed simultaneously. Simultaneously until first success: Simultaneous execution, but as soon as one action succeeds (i.e. has been completed or the phone is picked up), all others are terminated. Consecutively: All actions are executed in the specified order. Consecutively until first success: Consecutive execution, but as soon as one action succeeds, the following actions are not executed. Consecutively until first failure: Consecutive execution, but as soon as one action fails, the following actions are not executed.
Actions	Value	Explanation
Action 1	IP Notify: IQM VM2	Action Type and Profile: Select the Action Profile to be executed.
<input type="checkbox"/> Delete	0	Action Timeout or Duration: If this action runs longer than the time specified (0..3600 s), it is aborted and returns an error; 0 to deactivate. For Image Profile action, this is the duration and no error returns.

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MOBOTIX Community
Artikel zur IQ Messenger
Alarmintegration

NurseAssist Integrations Guide #4

Ascom SmartSense Nurse Call System

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The logo for Ascom, featuring the word "ascom" in a bold, lowercase, red sans-serif font.

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The screenshot shows the configuration page for 'IP Notify Profile 8' with the name 'SmartSenseFallDetected'. The page includes a 'Delete' button and a table with the following configuration details:

Profiles & Options	Value	Explanation
IP Notify Type	HTTP/HTTPS	Predefined Configuration: "MxCC Alarm" sends predefined network messages to the MxCC alarm list. <i>Acknowledge Required</i> prompts the MxCC user to confirm the message. If the alarm is not acknowledged within the specified acknowledge time, the camera triggers a transmission error. Select <i>Custom Configuration</i> to see the extended configuration.
Destination Address	dev-smartsense.ascom.com	Destination Addresses: Receiver IP address and port. Separate IP address and port using a colon. Enter one address per line.
HTTP Method	GET	HTTP Method: Transfer data using one of these HTTP methods.
	/api/thirdparty/alerts/trigger-alert?SensorId=c71DEYM&AlertNar	CGI-Path: Absolute CGI path beginning with '/'. This parameter allows using variables .
Security	HTTPS	HTTP Schema: Transfer data using these HTTP schemas.
	Basic authentication	Authentication Method: Method for HTTP Authentication.
	HTTP Authentication: User name and password for HTTP authentication separated by colon. Example: admin:meinSm
Data Type	JSON	Notification Data: Select type of IP notification data.

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```
/api/thirdparty/clients/status?SensorId={{SensorId}}&SensorTy  
atusName={{StatusName}}
```

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MOBOTIX Community Artikel
zur Ascom SmartSense
Integration

NurseAssist Integrations Guide #5

HPS ConectedHealth Nurse Call System

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The screenshot shows the configuration page for 'IP Notify Profile 9' with the name 'ConnectedHealth'. The interface is divided into two columns: 'Profiles & Options' and 'Explanation'. The 'Profiles & Options' column contains several fields with their current values:

- IP Notify Type:** HTTP/HTTPS
- Destination Address:** api.demo.connected-health.nl:3000
- Send Order:** Parallel send to all
- HTTP Method:** POST
- CGI-Path:** /kepler/event
- Security:** HTTP/1.1
- Authentication Method:** Basic authentication
- HTTP Authentication:** [Redacted]
- Data Type:** JSON
- Send Port:** 0

The 'Explanation' column provides detailed instructions for each field, such as 'Predefined Configuration', 'Destination Addresses', 'Send Order', 'HTTP Method', 'CGI-Path', 'HTTP Schema', 'Authentication Method', 'HTTP Authentication', 'Notification Data', and 'Port Number'. The 'Message' field is pre-filled with a JSON structure using variables like `$(MSG.LOCAL.KeplerNurseAssist.notification.uuid)`.

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General Settings	Value	Explanation
Action Group	FallDetected	Name: The name is purely informational.
	Enabled	Arming: Controls this action group: Enabled: activate the group. Off: deactivate the group. SI: group armed by signal input. CS: group armed by custom signal as defined in General Event Settings .
	(No time table)	Time Table: Time table for this action profile (Time Tables).
Event Selection	<input type="checkbox"/> Image Analysis: A5 <input type="checkbox"/> Image Analysis: VM <input type="checkbox"/> Image Analysis: VM2 <input type="checkbox"/> Message: NotInBed <input checked="" type="checkbox"/> Message: FallDetected	Event Selection: Select the events which will trigger the actions below. Use [Ctrl]-Click to select more than one event. Events in parentheses need to be activated first.
Action Details	5	Action Deadtime: Time to wait (0..3600 s) before a new action can take place.
	Simultaneously	Action Chaining: Choose how the status of each subaction influences the execution of all others. Simultaneously: All actions are executed simultaneously. Simultaneously until first success: Simultaneous execution, but as soon as one action succeeds (i.e. has been completed or the phone is picked up), all others are terminated. Consecutively: All actions are executed in the specified order. Consecutively until first success: Consecutive execution, but as soon as one action succeeds, the following actions are not executed. Consecutively until first failure: Consecutive execution, but as soon as one action fails, the following actions are not executed.
Actions	Value	Explanation
Action 1	IP Notify: ConnectedHealth	Action Type and Profile: Select the Action Profile to be executed.
<input type="checkbox"/> Delete	0	Action Timeout or Duration: If this action runs longer than the time specified (0..3600 s), it is aborted and returns an error; 0 to deactivate. For Image Profile action, this is the duration and no error returns.

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MOBOTIX Community Artikel
zur HPS ConnectedHealth
Integration

NurseAssist Integrations Guide #6

Skyresponse Alarm Nurse Call System

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MOBOTIX c71 Mobotix-c71-KVEC IP Notify Profiles

IP Notify Profile 10 SkyresponseFallDetected Delete

Profiles & Options	Value	Explanation
IP Notify Type	HTTP/HTTPS	Predefined Configuration: "MxCC Alarm" sends predefined network messages to the MxCC alarm list. Acknowledge Required prompts the MxCC user to confirm the message. If the alarm is not acknowledged within the specified acknowledge time, the camera triggers a transmission error. Select Custom Configuration to see the extended configuration.
Destination Address	sandbox.skyresponse.com	Destination Addresses: Receiver IP address and port. Separate IP address and port using a colon. Enter one address per line.
	Parallel send to all	Send Order: Send notification to one or more destinations. Sequential and parallel will send a notification to each destination address. Send to next on error will stop after the first successful notification or will try the next address if unsuccessful.
HTTP Method	POST	HTTP Method: Transfer data using one of these HTTP methods.
	/api/alarms	CGI Path: Absolute CGI path beginning with '/'. This parameter allows using variables.
Security	HTTPS	HTTP Schema: Transfer data using these HTTP schemas.
	None	Authentication Method: Method for HTTP Authentication.
		HTTP Authentication: User name and password for HTTP authentication separated by colon. Example: admin:meinSm
Data Type	JSON	Notification Data: Select type of IP notification data.
	<pre>{ "authentication": { "identifier": "XXXXXXXXXX", "secret": "39yAaWeyfwPJWMjg", "product": "KPLR" } }</pre>	Message: Message to include in Plain text notification data. When using HTTP protocol this text is used for QUERY_STRING in GET request. This parameter allows using variables.
Send Port	0	Port Number: Send a message from this camera port (0 for automatic).

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Art	Skyresponse-Kennung
CAMERA_UNREACHABLE_DETECTION	77431
HALLWAY_WANDER_DETECTION	77432
IN_BATHROOM_DETECTION	77425
INTRUDER_IN_ROOM_DETECTION	77433
MAN_DOWN_DETECTION	77421
MISSING_BED_DETECTION	77426
OUT_OF_BED_DETECTION	77423
OUT_OF_ROOM_DETECTION	77424
PERSON_VISIBLE_DETECTION	77427
PERSON_GETTING_UP_FROM_CHAIR_DETECTION	77428
PERSON_GETTING_UP_IN_BED_DETECTION	77429
SOEB_DETECTION	77422
STAFF_ENTERING_ROOM_DETECTION	77430
NO_DETECTION	77434
IN_BED_DETECTION	77434

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MOBOTIX AG
Kaiserstrasse
67722 Langmeil
Germany

+49 6302 9816-0
info@robotix.com
www.robotix.com

MOBOTIX, the MOBOTIX Logo, MxControlCenter, MxEasy, MxPEG, MxDisplay and MxActivitySensor are trademarks of MOBOTIX AG registered in the European Union, the U.S.A. and in other countries • Subject to change without notice • MOBOTIX do not assume any liability for technical or editorial errors or omissions contained herein • All rights reserved • © MOBOTIX AG